

PRODUCT SPECIFICATION

A palm-sized, easy-to-grip pendant for patients to activate calls on a TACERA Nurse Call System. The casing is manufactured from white, durable, impact-resistant plastic.

The 20 mm diameter button is easy to press, coloured green (or orange or red) for quick identification and contains dots for the letter "N" for the sight impaired.

An optional air pad pendant is available, which is a large 80 mm diameter soft touch, air-activated pad with cable for use by disabled patients or patients with reduced manual dexterity.

When a call has been initiated from the CALL button on the pendant, the wall callpoint beeps and a solid-state reassurance light within the CANCEL button on the wall callpoint activates. These combine to reassure your patient that the pendant is working correctly and that a call has been activated.

The callpoint CALL button and pendant removal alarm can be configured as different and independent calls.

The pendant cable is 2 m in length and manufactured from 3-core multi-strand, double insulated, white PVC cable and includes a linen clip that helps keep the pendant within easy reach of your patient. Longer cable lengths are available as an option.

The cable has a strong, moulded plastic, right-angle plug that connects to any TACERA callpoint with a 6.5 mm pendant socket. The design protects the callpoint socket, pendant plug and cable from damage should the cable be placed under severe tension, such as moving the bed with the pendant still connected.

The pendant electronics are protected against power surges, near lightning strikes, cable shorts and power reversals.

CISPR 22 Class B certification ensures that the pendant meet the new international emission requirements for hospitals and aged care facilities.

The pendant complies with the international RoHS Directive for the restriction of hazardous substances in electronic equipment.

Dimensions & Specifications

Height:	95 mm	Depth:	N/A
Width:	40 mm diameter	Weight:	165 grams
Housing:	Plastic moulding	Mounting Method:	Connects to any Austco callpoint with an optional 6.5 mm pendant socket
Connectors:	6.5 mm stereo jack	Cable Size:	2 m of 3-core multi-strand, double insulated, white PVC cable with a plastic linen clip attached
Operating Temperature:	0 – 50°C	Relative Humidity Range:	0 – 85% Non Condensing

Ordering Options

Code	Description	Order Options
CM-PM/xx	Pendant, Momentary Action, 2 m	xx: Specify button colour: GB – Green Button OB – Orange Button RB – Red Button Cable Length: For extra cable length, add /3M for 3 m /4M for 4 m /5M for 5 m /6M for 6 m /7M for 7 m e.g. CM-PM/xx/4M
CS-PMPAD	Pendant with Air Pad	

Designed to comply with:

- AS 3811: Hard-wired patient alarm systems.
- HTM 08-03: Bedhead services.
- UL 1069: Hospital Signaling and Nurse Call Equipment.
- IEC 60601-1, UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90: Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.



Manufactured in general accordance with the requirements of international quality assurance standard ISO9002.

Authorised Austco Reseller

