

# STERIBUTTON CALLPOINTS



Flush wall-mounted callpoints that initiate calls on the TACERA Nurse Call System. Back-lit buttons for night viewing by staff, reassurance LED and tone sounder for audibly or visually impaired patients, and call mute mode for cleaners are standard. Callpoints are also available with a 6.5 mm pendant socket or an 8 pin DIN socket for a patient entertainment handset.

*IPnet compatible and net- or web-configured*

*Flush soft touch silicone buttons with anti-bacterial additive*

*Easy to activate*

## LOWER "WHOLE-OF-LIFE" COSTS

The SteriButton series of callpoints are IP-based and feature plug-and-play functionality. Callpoints are connected in series using only standard CAT5/6 cable that also provides POIPnet (Power Over IPnet).

Each callpoint has a unique MAC address and is supplied configured for immediate use. Once connected, the callpoint is quickly and easily assigned to a room using a web-browser, which simplifies the process and reduces the cost of installation and commissioning.

The system continuously checks the functionality of each callpoint and will immediately notify maintenance staff if any problem is detected. SteriButton callpoints are hot-swappable, making maintenance a breeze – simply replace a callpoint and then use a web-browser to quickly re-assign functionality to the new device.

## EASY TO IDENTIFY AND USE

SteriButton callpoints feature soft-touch push buttons that are clearly labelled with text or icons for easy and fast identification. The call buttons are available in a number of standard colours and icons, allowing the system to operate

a variety of alarm levels and system inputs as determined by the protocols of your facility.

All callpoints are available with an optional socket for a Call Pendant or Patient Entertainment Handset.

Users can easily initiate a call by pressing the call button on the callpoint or from an optional plug-in Call Pendant.

All callpoints are connected on the IPnet patient area network to all other components of the TACERA system to initiate communication between nursing staff and patients.

## IMPROVED INFECTION CONTROL

These ergonomic callpoints are flush, wall-mounted with a configurable CALL button. Each callpoint has a separate CANCEL button that illuminates red, giving visual reassurance that a call has been initiated. Both buttons are manufactured from silicone rubber with an anti-bacterial additive for improved infection control.

To avoid nuisance alarms during cleaning, the callpoint incorporates a cleaning mode. The cleaning time and date is logged on the central server, providing an important audit trail for infection control compliance.

- A pendant removal alarm ensures instant notification that a Call Pendant is not available for your patient, reducing downtime.
- SteriButton callpoints fit into standard electrical wall boxes, reducing installation times and costs.
- Manufactured from an impact-resistant polycarbonate / ABS plastic blend, all callpoints feature concealed fixings to maximise the life of the callpoint and reduce maintenance costs.

- For vital Emergency and Code Blue alarms, the callpoint can be configured to show elapsed time on the annunciator and other devices.
- Gold plated button-contacts prevent corrosion and ensure the callpoint is functional, maximising the life of the callpoint and reducing maintenance costs.
- All callpoints buttons are AS 3811 compliant.
- Class B certification ensures the callpoints meet the new emission requirements for hospitals and aged care facilities.